



Product information | TECHNICAL DATA SHEET

ETHYL ALCOHOL 70%

400 ml, 10.86 Oz, 308 g. 100 ml, 2.64 Oz, 75 g.

GENERAL DESCRIPTION

Hydroalcoholic spray solution ideal for cleaning surfaces of all kinds.

FEATURES

- Degreasing power.
- Ease of application.
- Hydroalcoholic Solution.
- Multi-surface cleaner. No added perfumes.

APPLICATIONS

Hydroalcoholic solution ideal for cleaning surfaces of all kinds, including glass, metals and plastics. Ideal for methacrylate: helmet visors, protective covers, light screens, etc.

In vehicles: steering wheels and other interior and exterior plastics. Home furniture and all office items, such as tables, screens and computer keyboards or mobile phones. The product may be used to easily remove residues from grease, ink stains, paint and traces of glue from self-adhesive labels from all types of surfaces. Due to its formulation, it may provoke reactions on some inks and paints, therefore it is recommended to do a preliminary test to verify that the product does not compromise the aesthetic appearance of your objects.

INSTRUCTIONS FOR USE

- Shaking is not necessary.
- Use the can in a vertical position and apply at a distance of about 20 cm.
- The surface is to be cleaned must be dampened with the product and left to evaporate. On screens or glass, allow for the product to act for a minute and remove it with a soft cotton cloth or blotting paper.
- Do not apply to electric equipment that is plugged in.
- Request a safety data sheet.

PHYSICAL PROPERTIES

Alcohol/water ratio:	70% (v/v)
Valve:	Femal
Theoretical yield (continuous painting):	15 m ² for 400 ml
	3,75 m ² for 100 ml

Montana Colors, S.L. P.I. Clot del Tufau - 08295 Sant Vicenç de Castellet BARCELONA

T.: +34 938 332 787 - montana@montanacolors.com - www.montanacolors.com









APPEARANCE



LID

Protects the valve and the nozzle.

LABEL CAP SYSTEM





Pocket Cap

EX014PR0003 Ethyl alcohol 70%





PACKAGING

Nominal capacity: 100 - 400 ml. Heat resistance of packaging: <50°C Propellant: DME

Box (6 units) - 400 ml 205 X 140 X 200 mm.

Box (24 units) - 100 ml 275 X 190 X 150 mm.









1/1